
Z21047104

User Manual

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CryalJet® Maxi/Mini

Distribution:
AIR LIQUIDE Parc Gustave Eiffel, 8 avenue Gutenberg,
Bussy Saint Georges, 77607 Marne la vallée Cedex 3, France

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FILLING THE UNIT

The CryoPro® unit is to be filled with liquid nitrogen only. Fill the unit by slowly pouring liquid nitrogen into the top or by using a standard low pressure withdrawal device from a storage tank (dewar). Depending on the number of lesions to be treated and the duration of the individual treatments the unit should be filled to 25 - 90 % of its maximum capacity. We recommend the unit to be approx. 70 % filled for a 4 - 6 hour duration of intermittent use. After filling a warm unit, the liquid will boil violently for approx. 30 secs. Hereafter, the top cover can easily be replaced. A weak but audible hissing sound from the pressure relief valve is normal and does not indicate an error condition. Caution should be taken when refilling a CryoPro® containing residual fluid as the unit is pressurized at 0.7 bar (10 psi). Slowly unscrewing the top cover will allow the unit to depressurize before the top comes off. This is accompanied by an audible hissing sound until it is fully depressurized.

OPERATION

The CryoPro® is supplied with a set of standard spray apertures. Optional closed probes are available for contact freezing. All probes and spray apertures are screwed into the permanently affixed knurled nut with fingertight firmness. Do not apply excessive force.

Do never use the CryoPro® without either a spray aperture or closed probe in place. The excessive amount of liquid put out if the unit is used without spray or probe attachment will cause severe damage to the patient.

Spray freezing (REF OS A, B, C, D, SS, BS, ACNE, LL)

The standard spray apertures consist of: "A" (1 mm opening) and "B" (0.75 mm opening) for large and deep lesions or for superficial desquamation of skin; 2 X "C" (0.5 mm opening) recommended for warts and keratoses; "D" (0.37 mm opening) for small warts and keratolytic lesions. Finally, the "Bent Spray" aperture (0.5 mm opening) is for areas difficult to reach like behind the ears, lip and outer nose (avoiding gas into the airways). Select a spray aperture appropriate for the lesion to be treated. In-depth freezing is best obtained when the spray aperture is close to the lesion (5 - 10 mm distance) as the liquid content of the spray is higher close to the aperture. When possible, raised warts should be approached tangentially from at least two sides. This allows the ice ball to move down through the wart while sparing the surrounding tissue. Spraying directly at the lesion (perpendicular to the surface) will cause more lateral spread with less penetration. To avoid this apply the spray intermittently or use a smaller spray aperture. Superficial desquamation is obtained using the larger apertures ("A" or "B") at approx. 5 cm distance (higher vapor content) and slowly "painting" the surface with vaporized nitrogen.

Contact freezing (REF CP SP, CP Imm-30mm, CX DS19, CX DS25, CX 19x5, CX 19x25, CX 25x5)

Closed contact probes provide in-depth freezing with minimal lateral spread. Select a probe size fitting the size of the lesion. Make the lesion wet with a drop of water or contact gel prior to freezing to facilitate the thermal transfer between skin and probe tip. Apply the contact probe when still warm applying light pressure and pull the trigger handle. Allow approx. 5 sec. extra freezing time to cool down the contact probe and keep freezing until a 1 mm halo is seen around the probe tip. Then stop the freezing and allow a few seconds for the probe tip to come loose. This "quick freeze slow thaw" procedure using the contact probe technique ensures a very effective treatment with high penetration and accuracy.

Freezing times vary according to the aperture or probe used, the size and type of the lesion and the freezing distance (spray only). The user is encouraged to familiarize himself with available literature on cryosurgery using liquid nitrogen.

AFTER USE

At the end of an office day it is highly recommended to store the CryoPro® with the lid on in order to prevent condensation inside the bottle. Build-up of condensation may lead to clogging of the valve assembly.

MAINTENANCE

No preventive maintenance is needed.

STERILIZATION

Normally, spray apertures need no sterilization as they are not in contact with the patient. However, should the need for sterilization arise, the same procedure as for closed contact probes can be followed. Closed contact probes are in direct contact with the patient and should always be sterilized after use. The recommended method for sterilization is steam autoclaving at 121 deg. C for 15 min.

WARNINGS

- Do never use the CryoPro® unit without a spray aperture or closed contact probe in place.
- The CryoPro® is supposed to be held in a vertical position when used.
- Liquid nitrogen may escape through the relief valve if the unit is held in a position allowing the liquid inside the bottle to be in contact with the lid. Return the unit to a more upright position if liquid starts coming out of the relief valve. The maximum working angle depends on the amount of liquid nitrogen in the bottle.
- Be careful when opening a pressurized unit. Unscrew the top slowly. The unit may still be pressurized even if there is no residual liquid left.
- The CryoPro® is a surgical instrument. When not in use it must be kept out of reach of children and unauthorized personnel.

STORING LIQUID NITROGEN

To ensure a clean supply of liquid nitrogen the storage tank (dewar) must be completely emptied prior to refilling 3 - 4 times per year. If necessary, water from ice crystals remaining inside the tank can be dried out by creating air circulation (e.g. using a vacuum cleaner for approx. 15 min.).

SUGGESTED LITERATURE

- Dawber, Colver, Jackson: Cutaneous Cryosurgery, Principles and Clinical Practice (Martin Dunitz ISBN 1 85317 082 8).
- Kuflik, Gage: Cryosurgical Treatment For Skin Cancer (Igaku-Shoin ISBN 0 89640 157 X)

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CONGELAMENTO SPRAY (REF. OS A, B, C, D, SS, BS, ACNE, LL)

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CONGELAMENTO A CONTATTO (REF. CP SP, CP 1mm-30mm, CX DS19, CX DS25, CX 19x5, CX 19x25, CX 25x5)

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LETTERATURA CONSIGLIATA

- Dawber, Colver, Jackson: Cutaneous Cryosurgery, Principles and Clinical Practice (Martin Dunitz ISBN 1 85317 082 8).
- Kuflik, Gage: Cryosurgical Treatment For Skin Cancer (Igakusho-in ISBN 0 89640 157 X)

LN₂ Treatment Times	Technique	Time (s)	FTC	Margin (mm)	Sessions/interval (weeks)
Actinic keratosis	OS	5	1	-	- / -
Basal cell carcinoma (1-1.5 cm)	OS or CP	40-90	1	3-5	1 / -
Cervical erosions	CP	45-60	1	-	1 / -
Cherry angioma	CP	10	1	-	1 / -
Condylomata accuminata	OS or CP	5-10	1	-	2-4 / 3-4
Dermatofibroma	CP	20	1	1-2	1-2 / 8-10
Lentigo, solar + simplex	OS	5	1	-	1 / -
Molluscum contagiosum	OS or CP	5	1	-	1 / -
Nevi, epidermal	OS or CP	5	1	1	2 / 6
Seborrhoeic keratosis	OS	5-10	1	1-2	1 / -
Skin tags	OS	5	1	1	1 / -
Verrucae planae	OS	5	1	1	1-2 / 3-4
Verrucae seborrhoicae	OS or CP	5-10	1	-	1 / -
Verrucae vulgaris	OS	10-15	1	1-2	1-3 / 3-4
Veruccae plantaris	OS	10-20	1-2	1-2	2-4 / 3-4

OS = open spray CP = contact probe FTC = freeze/thaw cycles

Treatment times are indicative and for reference only and may vary according to the size of the spray aperture or closed contact probe used, the size and type of the lesion, the freezing distance (spray only) as well as personal preferences. The user is encouraged to familiarize himself with spray freezing and contact probe freezing techniques as described in the user manual. ...CryoPro® & CryAIJet® by CORTEX TECHNOLOGY

Customer satisfaction survey

It is our goal to ensure the highest possible satisfaction rate by providing quality products that fulfill the needs of our customers. Accordingly, you are encouraged to complete the feedback form below and return your comments by fax or to go to our feedback web-page <http://www.cortex.dk/feedback.htm>. Your assistance is highly appreciated.

Product name:

Model:

Config. (if applic.):

Comment on performance, function, design and/or ease of use:

Report problems and malfunction below:

Other issues:

Name:

Email:

Fax to: +45 9857 2223 or visit <http://www.cortex.dk/feedback.htm>

0543